

cs@american-sp.com +1-734-884-4770american-sp.com 17177 N. Laurel Park Dr. Suite 340 Livonia, MI 48152

Urinalysis Tubes

RODUCTS

AMERICAN

The American Scientific Products Urinalysis Tube is intended for use in the protection of fresh urine specimens from the time of voiding, through holding and transport to the clinical laboratory in the temperature range from 4 - 25 degrees C. for a time period not to exceed 72 hours. The ASP Urinalysis Tube preservation system allows room temperature exposure of the urine specimen (Volume: 2-0ml) for up to 3 days without alteration of reagent strip chemistries or microscopic examination and/or enumeration of formed elements. Made with copolymer containing built-in antimicrobial preservative. Tube comes with cap (separate).



Cat. No.	Description	Unit
ASP-8598-QUS1P	Tube UA Trans White Cap No Presv., Printed	1000/Case
ASP-9000QW	Tube UA W/ White Cap & Persv. Blue Lines	1000/Case
ASP-9000Q	Tube UA W/Cap & Persv. Blue Lines	1000/Case
ASP-9100Y	Tube, Urinalysis Preserv 10ml w/Yellow top Attached	1000/Case
ASP-8589TUY	Caps for urine tubes yellow heavy	1000/Case

Ordering Information

